

User's hints

BGE constituents and analytes hints:

- Sort constituents by name or by any of their displayed parameters via clicking on the relevant column header of the *BGE Constituents* or *Analytes* subwindow.

Electropherogram hints:

- Point to the individual peaks with the mouse cursor to read the names of the corresponding analytes in a hint window.
- Quick-zoom your chart via drawing a rectangle (start at the upper left corner) around the area you want to see in detail.
- Undo your quick-zoom via double clicking on the chart.

Database of constituents hints:

- Click *Pick from database...* button in the *Add BGE constituent...* or *Add analyte...* window (appears after clicking the *Add...* button on the main form) to retrieve values of a database constituent.
- Add, remove, or edit database constituents via clicking *Tools|Modify database...* from the main menu. However, modifying the values of database constituents is recommended only if your modifications have general validity, i.e. if they are true always, not only in the current document or configuration. If you need modified constituents only for the current document or configuration, you are invited to add them in the *BGE Constituents* or *Analytes* subwindow as a new constituent (under a new name, e.g. "sodium-modified").